





MicroBlast™



***Automated Deionizing
Air Knife Plate Cleaner***

MicroBlast™

Automated Deionizing Air Knife Plate Cleaner



MicroBlast is an integratable high-throughput device for cleaning labware with a deionizing blast of air. MicroBlast can even clean lidded plates; the instrument will automatically lift the lid and clean separately, then relid the plate. MicroBlast can be used stand-alone or with an automated system.

MicroBlast Key Features:

- Cleans top and bottom of any microplate
- Cleans the lid of lidded plates
- Compatible with MicroCart™ and MicroDock™
- Automation friendly with easy-to-use integrators' kit
- Stackable

MicroBlast Key Specifications:

Dimensions (H x W x L)	12.6" x 14.2" x 27.2"
Compressed Air	@80psi, 0.004cfm
Power	120/240VAC, 50/60Hz
Weight	88 lbs

The Building Blocks for Discovery.

HighRes Biosolutions

299 Washington Street

Woburn, MA 01801

Tel 781.932.1912

Fax 781.938.0813

www.highresbio.com