



MicroPierce™



Plate Piercer

MicroPierce™

Plate Piercer



MicroPierce is the smallest and most functional plate piercing instrument available. MicroPierce features online pin washing and multiple pierce heights, and can pierce deep or shallow well plates.

MicroPierce can pierce up to 20 foil seal layers, completely opening each well for non-contact dispensing or pipetting. Compact for integration into an efficient laboratory setup, the MicroPierce tool fits well into your robotic process and offers the added convenience of a pin-washing feature.

Key Features:

- Nearly half the size of its closest competitor
- Features online pin washing
- Ability to pierce multiple foil layers
- 384 and custom pierce heads available
- Compatible with MicroCart™ and MicroDock™
- Automation friendly with easy-to-use integrators' kit
- Stackable

Dimensions (H x W x L)	8.5" x 8.9" x 12"
Compressed Air	@80psi, 0.06cfm
Power	120/240VAC, 50/60Hz
Weight	29 lbs

The Building Blocks for Discovery.

HighRes Biosolutions

299 Washington Street

Woburn, MA 01801

Tel 781.932.1912

Fax 781.938.0813

www.highresbio.com