




NanoServe



High-Density Labware Handling



NanoServe™

High-Density Labware Handling

NanoServe is a high-density way to deliver plates and tips to your automated systems, handling up to 220 standard microtiter plates in a small footprint. NanoServe is efficient and compact for high-throughput applications. Plates stack one on top of the other in removable stackers, so loading and unloading are easy. The stackers are compatible with MicroServe. NanoServe works with most SBS footprint labware.

Key Features:

- Load 220 standard microplates in an incredibly small footprint
- Easy to load and remove stackers
- Compatible with MicroCart™ and MicroDock™
- Automation friendly with easy-to-use integrators' kit

Dimensions (H x W x L)	36" x 25" x 21"
Power	120/240VAC, 50/60Hz
Weight	100 lbs (empty)
Capacity	700mm available for stacking per hotel, 4 hotels Up to 220 standard height microtiter plates Up to 44 tip boxes

The Building Blocks for Discovery.

HighRes Biosolutions

299 Washington Street

Woburn, MA 01801

Tel 781.932.1912

Fax 781.938.0813

www.highresbio.com